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Product description/Use

The main function of the plastic hose clamp is to close the connection tubes of medical equipment, creating a physical barrier for the liquids contained in the tubes. The closing of the tube is assured and maintained by a catch located at the contact point of the two gripping rings. The catch consists of two plastic strips both with notches to be able to clasp the device in different positions. On the basis of the notch position selected, less or more pressure will be exercised on the tube or the material collected. The rear part, between the gripping rings and the central pin, has a knurled area which may be used as flow limiter by inserting the tube and clasping the clamp.

The device may be used on PVC tubes having a hardness ranging between 60 and 88 ShA and on silicone tubes.

The hose clamp is made of Nylon reinforced with 30% glass fiber. Latex free. Phthalate free.

Sterilizable with: Ethylene Oxide; Gamma (up to 25 kGy); Steam (cycle at 125°C for 20 minutes).

Warnings and precautions

- ❑ For delicate tubes (ex. silicone tubes), do not fully tighten the device, but clasp it in the first three positions of the rear catch.
- ❑ After having fixed the device, exercise a slight outward opening force to check the tightness of the catch.
- ❑ Place the tube in central position taking care to never insert it right to the bottom in contact with the fulcrum to prevent excessive deformation of the tube.
- ❑ Always wear the protective gloves to prevent any type of biological risk during the use of the device and risks from contact with biological fluids or with substances being used in the tubes.
- ❑ Always dispose after use.
- ❑ If the use is different from that suitable, it is the buyer the only responsible for the choice, the adequacy and the use of this product.
- ❑ **It is necessary that the user and/or patient notify any serious incident connected to the device, manufacturer and competent authority of the state in which the user and/or patient is established.**

Side effects

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SA.GE Srl - Via Montale 110 – 41032 Cavezzo (MO) Italia

info@sage-srl.com – www.sage-srl.com

No particular side effects resulting from the use of the device have been found.

No interferences or risks, if used along with other devices, have been found.

Disposal

The device must be disposed of as "Hospital waste potentially contaminated".

THE RISKS FOR THE PATIENT RESULTING FROM THE REUSE OF THE DEVICE: THERE ARE THOSE OF CROSS-CONTAMINATIONS AND INFECTIONS, INFECTIVE COMPLICATIONS AS LOCAL INFECTIONS AND SYSTEMIC INFECTIONS.

Storage and Handling

Store the device in a dry and clean place away from direct heat source.

Replace the device if it does not appear completely intact.

NOTE: It is essential to read these instructions for proper use of the device.

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