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**EN**

**CODE AC 165 FISTULA ARM CLAMP**  
Class I non sterile medical device, without measuring function.

**CE****Intended use**

Device to be used as an aid for dressing procedures (injections, blood sampling, hemodialysis, etc ...)

**Cautions**

Check the integrity of the device before use.

The device must be used by personnel qualified to perform generic medication operations.

Apply the personal safety rules when handling potentially infected devices after use on the patient.

Do not use devices found unattended, dirty, used and not disinfected, etc ...

**Instruction for use**

- Properly held the device resting the external handle on the palm of the hand and gripping the internal handle with the other fingers above the spring.
- Apply slight force on the grip to open the device.
- Insert it on the site involved adequately medicated.
- Slowly release and close the device making sure that compression is applied on the site involved in the injection or blood sampling.

**Cleaning and disinfection**

The device is reusable. After use, the device must always be washed and disinfected. First wash and then disinfect the device by immersing it in enzymatic detergent and then disinfecting solution based on peracetic acid (for example Viroxid). Follow the instructions/dosages of the detergent-disinfectant used. Store in a sealed bag.

**Warning**

- Before using the device, carefully medicate the site involved in the injection or sample-taking disinfecting and placing cotton wool and gauze between the epidermis and the device.
- Never place the device in direct contact with the wound.
- Do not use the device to compress cuts, lacerations, etc. or in any case for purposes other than the intended use.
- Wear protective gloves in order to prevent any biological risk during use of the device and risks due to contact with biological fluids or other substances that could be dangerous.
- Dispose the device if contaminated with physiological fluids, blood and/or other agents.
- It is necessary that the user and/or patient notify any serious incident connected to the device to the manufacturer and the competent authority of the state in which the user and/or patient is established.**

**Contraindications and adverse reactions**

Use of the device does not produce any side effects.

No interferences or risks in use associated with other devices have been reported.

**Discard after use**

Dispose the device as medical hospital waste potentially contaminated.

**Storage conditions and handling**

Store the device in a dry and clean place away from direct heat sources.

Replace the device if it is not fully intact or if its essential parts show signs of wear or if the spring has lost its functionality.

**NOTE:** it is essential to read these instructions for proper use of the device.

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**SA.GE Srl** - Via Montale 110 – 41032 Cavezzo (MO) Italia  
[info@sage-srl.com](mailto:info@sage-srl.com) – [www.sage-srl.com](http://www.sage-srl.com)

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