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Class I non-sterile medical device.

**Product description/Use**

The Multifunction Clamp can be used for various medical uses.

Its main use is to be a means of grasping gauze and tampons during generic dressing operations. It can also be used to stop sheets.

Clamp closure is ensured and maintained by a hook placed at the point of contact of the two rings for the handle. The attachment consists of two plastic flaps both with recesses to be able to hook the device in different positions. According to the chosen position, a different pressure is exerted on the pipe or on the collected material, lower in the first hooks and higher in the last positions.

Made of Nylon with 30% glass fiber. Latex free. Phthalate free.

Sterilizable with: Ethylene Oxide; Range (up to 25kGy); Steam (cycle at 125° for 20 minutes).

Warnings and precautions

- ❑ After having fixed the device, exercise a slight outward opening force to check the tightness of the catch.
- ❑ During the use of the device, always wear the protective gloves to prevent any type of biological risk.
- ❑ Always dispose after use.
- ❑ If the use is different from that suitable, it is the buyer the only responsible for the choice, the adequacy and the use of this product.
- ❑ **It is necessary that the user and/or patient notify any serious incident connected to the device to the manufacturer and the competent authority of the state in which the user and/or patient is established.**

Side effects

No particular side effects resulting from the use of the device have been found.

No interferences or risks, if used along with other devices, have been found.

Disposal

The device must be disposed of as "Hospital waste potentially contaminated".

THE RISKS FOR THE PATIENT RESULTING FROM THE REUSE OF THE DEVICE: THERE ARE THOSE OF CROSS-CONTAMINATIONS AND INFECTIONS, INFECTIVE COMPLICATIONS AS LOCAL INFECTIONS AND SYSTEMIC INFECTIONS.

Storage and handling

Store the device in a dry and clean place away from direct heat source.

Replace the device if it does not appear completely intact.

NOTE: It is essential to read these instructions for proper use of the device.

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