**IT****CODICE AC 170 PINZA ANATOMICA**

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Class I non-sterile medical device.

**Product description/Use**

Device to be used to grasp gauze and tampons during general dressing operations.

Latex free

Phthalate free

**Usage precautions**

Check the integrity of the device in all its parts before using it.

The device must be used by qualified personnel to intervene for generic dressing operations.

Apply the personal safety policy when handling potentially infected, patient-originating devices.

Absolutely avoid the use of devices removed or found unattended, dirty, used, etc.

**Warnings and precautions**

- The device must not be used on injured skin, it must not be used directly in contact with the skin.
- Discard after use, especially if it is contaminated with physiological fluids, blood and / or other agents.
- Replace the device after each use.
- Do not sterilize the device in a steam autoclave.
- If the use is different from that suitable, it is the buyer the only responsible for the choice, the adequacy and the use of this product.
- It is necessary that the user and/or patient notify any serious incident connected to the device to the manufacturer and the competent authority of the state in which the user and/or patient is established.**

**Side effects**

No particular side effects resulting from the use of the device have been found.

No interferences or risks, if used along with other devices, have been found.

**Disposal**

The device must be disposed of as "Hospital waste potentially contaminated".

**THE RISKS FOR THE PATIENT RESULTING FROM THE REUSE OF THE DEVICE: THERE ARE THOSE OF CROSS-CONTAMINATIONS AND INFECTIONS, INFECTIVE COMPLICATIONS AS LOCAL INFECTIONS AND SYSTEMIC INFECTIONS.**

**Storage and handling**

Store the device in a dry and clean place away from direct heat source.

Replace the device if it does not appear completely intact.

**NOTE:** It is essential to read these instructions for proper use of the device.

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