

DATA SHEET

Code


Description

AC 150 B
AC 150 R
AC 150 V

Scissor tubing clamp – blue red green



Picture



Class I non-sterile medical device certified according to EU Regulation 2017/745 (MDR) according to annex II+III.

Manufactured in 30% glass fiber filled Nylon, to be used for the partial or total closure of plastic / silicone pipes with hardness between 60 and 80 ShA, 4 closing positions.

Latex free

Phthalate free

Quality control: visual

Sterilizable with: Ethylene Oxide; Range (up to 25kGy); Steam (125° cycle for 20 minutes)

Standard package: 800 pieces in bulk, not sterile, in double polyethylene bag and carton.

Conservation: Store the medical products in a dry and clean place, away from sources of direct heat.

Storage temperature: 0 ° C / 40 ° C (32 ° F / 104 ° F)

Self life: 7 years

Caution: Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

Registration: Code AC 150-CND V9099-GMDN 40589-N° Publication RDM 2206064