

DATA SHEET

Code



Description

**AC 155 B
AC 155 V**

Multipurpose scissor clamp – blue or green



Picture



Class I non-sterile medical device certified according to EU Regulation 2017/745 (MDR) according to annex II and III.

Manufactured in 30% glass fiber filled Nylon.

Device to be used to grasp gauze and tampons during general dressing operations.

It can also be used to stop sheeting. 5 closing positions.

Latex free

Phthalate free

Quality control: visual

Sterilizable with: Ethylene Oxide; Range (up to 25kGy); Steam (125° cycle for 20 minutes)

Standard package: 600 pieces in bulk, not sterile, in double polyethylene bag and carton.

Conservation: Store the medical products in a dry and clean place, away from sources of direct heat.

Storage temperature: 0 ° C / 40 ° C (32 ° F / 104 ° F)

Self life: 7 years

Caution: Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

Registration: Code AC 155-CND V9099-GMDN 47570-N° Publication RDM 2205984