


**DATA SHEET**

Code	Description	
AC 165 V AC 165 B AC 165 CO	Fistula arm clamp – green or blue	CE

**Picture**



**Class I non-sterile medical device certified in accordance with EU Regulation 2017/745 (MDR) according to annex II+III.**

Manufactured in 30% glass fiber filled Nylon, Aisi 320 stainless steel spring according to DIN 17224-1, load of 4.9 kg at maximum opening. To be used as an aid for dressing procedures (injections, blood collection, hemodialysis, etc ...).

Maximum opening of the clamp: 110 mm.

Latex free

Phthalate free

**Quality control:** visual

**Sterilizable** with: Ethylene Oxide; Range (up to 25kGy); Steam (125° cycle for 20 minutes)

**Standard package:** 50 / 100 pieces in bulk, not sterile, in double polyethylene bag and carton or carton of 25 neutral boxes containing 2 pieces.

**Conservation:** Store the medical products in a dry and clean place, away from sources of direct heat.

**Storage temperature:** 0 ° C / 40 ° C ( 32 ° F / 104 ° F )

**Self life:** 7 years

**Caution:** Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

**Registration:** Code AC 165-CND F9080-GMDN 446604-N° RDM Publication 2206043