


DATA SHEET

Code	Description	
AC 170 B (113 mm) AC 170 V (113 mm) AC170LV (138 mm)	Anatomic tweezer – green or blue	CE

Picture



Class I non-sterile medical device certified in accordance with EU Regulation 2017/745 (MDR) according to annex II+III.

Manufactured in polypropylene material with 5% glass fiber and 20% talc, to be used to grip gauze and tampons during general dressing operations.

Latex free

Phthalate free

Quality control: visual

Sterilizable with: Ethylene Oxide; Range – Beta (up to 25kGy)

Standard package: 1000 pieces in bulk, not sterile, in double polyethylene bag and carton.

Conservation: Store the medical products in a dry and clean place, away from sources of direct heat.

Storage temperature: 0 ° C / 40 ° C (32 ° F / 104 ° F)

Self life: 7 years

Caution: Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

Registration: Code AC 170-CND V9099-GMDN /-N° Publication RDM 2206079