


DATA SHEET

Code	Description	
AC175LV (138 mm)	Anatomic tweezer – green	CE

Picture



**Class I non-sterile medical device certified in accordance with EU Regulation 2017/745 (MDR) according to annex II+III.**

Manufactured in polypropylene material with 5% glass fiber and 20% talc, to be used to grip gauze and tampons during general dressing operations.

Latex free

Phthalate free

**Quality control:** visual

**Sterilizable** with: Ethylene Oxide; Range – Beta (up to 25kGy)

**Standard package:** 1000 pieces in bulk, not sterile, in double polyethylene bag and carton.

**Conservation:** Store the medical products in a dry and clean place, away from sources of direct heat.

**Storage temperature:** 0 ° C / 40 ° C ( 32 ° F / 104 ° F )

**Self life:** 7 years

**Caution:** Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

**Registration:** Code AC 175-CND V9099-GMDN-N° Publication RDM 405945