





DATA SHEET

Code	Description	
AC 180 B	Kelly type scissor clamp - blue	
CE		
Picture		
		
		

Class I non-sterile medical device certified in accordance with EU Regulation 2017/745 (MDR) according to annex II + III.

Device to be used to grasp gauze and tampons during general dressing operations.

4 closing positions.

Latex free

Phthalate free

Quality control: visual

Sterilizable with: Ethylene Oxide; Range (up to 25kGy); Steam (cycle at 125 ° for 20 minutes)

Standard package: 1000 pieces in bulk, not sterile, in double polyethylene bag and carton.

Conservation: Store the medical products in a dry and clean place, away from sources of direct heat.

Storage temperature: 0 ° C / 40 ° C (32 ° F / 104 ° F)

Self life: 7 years

Caution: Exposure of this article to materials, chemicals, sterilization modes or temperatures outside of specifications may cause changes in the product material.

If the use is different, the buyer is solely responsible for the choice, adequacy and use of this product.

Registration: Code AC 180-CND V9099-GMDN 47570-N° Publication RDM 2206071